

MEASUREMENT GUIDE FOR COLOPLAST FREEDOM® BRAND MALE EXTERNAL CATHETERS

Coloplast
Customer Service
888-726-7872

M1634N 12.07

with Kink Resistant Nozzle

Freedom Clear® LS Song Seal • Freedom Clear® SS Sport Sheath • Freedom Clear®

Freedom Cath® • Active Cath® • UroSan® Plus • Gizmo™
with Aloe and Kink Resistant Nozzle Clear Advantage®

X-Large-40 mm

Large-35 mm

Intermediate-31 mm

Medium-28 mm

Small-23 mm

Print this file and cut out the half-moon shapes. Measure yourself in a relaxed state, closest to your body. If you find yourself in between 2 sizes, manufacturer recommends you size down. See our FAQ's on any of the Male External Catheter pages for more information.